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Do not use when folded or creased

Do not insert pins

Not to be used by very young children (0-3 years)

Products:

Item No.: 93050 – FIR Infrared Knee Bandage

Item No.: 93051 – FIR Infrared Waist Bandage

Item No.: 93052 – FIR Infrared Wrist Bandage

System components and functions

The delivery comprises a rechargeable control battery (LBP-742D lithium-polymer battery pack), a FIR infrared bandage and a power adapter.

- The control element (rechargeable battery) is used to regulate the temperature.
- The FIR infrared bandage warms the body.

Getting started:

1. Charge the battery until the indicator (LED lamp) lights up in green. When the battery is empty, the LED lamp will light up in red. Note: The battery has to be fully charged before the first use.
2. First wipe the inner padding of your FIR infrared warmth bandage with a moist cloth to support the heat conduction.
3. Apply the bandage to the body region to be treated.
4. Place the battery in the dedicated bag on the bandage, and connect the cable of the bandage with the control element (battery).
5. Push the ON/OFF button. The LED light will light up in red to indicate that the device works on its maximum level of intensity.
6. Press the ON/OFF button again. The LED will light up in green to indicate that the device works on its minimum level of intensity.
7. Press the ON/OFF button once more to switch off the device.
8. 15-30 minutes after starting the device, you will feel a thermal reaction.
9. After using the device, wipe off the sweat with a moist cloth and let the bandage dry by placing it horizontally on a surface for at least one hour.
10. Carefully fold up the bandage and store it in a dry and cool place.

Caution!

1. This product is not a medical device and it may not be used for medical purposes.
2. We recommend to drink a glass of warm water before and after use to stimulate your metabolism.
3. In case you take remedies, ask your doctor whether you may use the bandage to be sure.
4. If your skin burns or turns persistently red after using the bandage, this will be a normal phenomenon. It will gradually decrease and finally disappear completely.
5. The length of application can be determined individually.

Warning!

1. Do not switch on the heating-up function when the bandage is folded up.
2. Keep the bandage from contact with sharp or pointed objects in order to prevent its outer material from cuts or other damage.
3. Do not use the bandage when the upper material has been damaged.
4. Do not use spare or additional components that are not expressly indicated as original parts.
5. Do not disassemble or repair the device or install spare parts into it as this will void the guarantee.
6. Handle the bandage with care whenever it is switched on, especially when there are children nearby.
7. Disconnect the power plug when the battery is fully charged.
8. Never use the bandage in water or in a moist environment.
9. Keep the bandage away from electrical sources of heat.
10. Do not throw or drop the bandage.
11. Keep the bandage out of the reach of children.
12. Persons who have difficulties using the bandage should seek support.
13. If the LED goes out during the use of the bandage, the battery needs to be recharged.
14. Disconnect the device from the power grid in case of unusual smells during the use or recharging process of the device.
15. Rechargeable battery: The battery should be charged every 90 days to ensure its functionality.
16. Warning: Please do not use the product if you are under the following conditions:
 - severe conditions of the heart or brain
 - when you are wearing implants or pacemakers, insuline pumps, etc.
 - local damage or injuries of the skin
 - during pregnancy and breastfeeding
 - in case of cancer or other critical diseases
 - allergies
 - acute traumas

Care and storage

1. Do not clean the bandage with benzene, diluters or other strong detergents.
2. Always use a dry cloth for cleaning the control battery – never use a moist or wet cloth for cleaning.
3. To ensure hygiene and maximum durability of the bandage, clean its surfaces with a moist cloth after the heating phase – a cotton cloth moistened with water would be best.
4. This bandage does not require any special maintenance.
5. If you do not use the bandage for a longer time, store it with its accessories in a dry and cool environment.
6. Never expose the product to high temperatures, e. g. never place it near fire. Never expose the product to direct sunlight for a longer time.

Protection of the environment

At the end of its life cycle, this product must not be disposed of in household waste due to its battery. In order to dispose of your device, take it to a collection point for the disposal of electric or electronic devices. This is what the symbol on the product, in the user's manual, or on the package hints at. The materials of the device can be recycled according to their labelling. By recycling material and returning discarded devices you are making a considerable contribution to environmental protection. Please ask your local authorities for the appropriate disposal point.

Technical data

Materials:

90 % polyester, 10 % polyamide

A. Charger model # LBC-8412

1. Input voltage: 100 ~ 240 V AC; 50/60 Hz
2. Charging voltage: 8.45-8.6 V DC
3. Charging current: 0.9-1.1 A

B. Rechargeable lithium-polymer battery: Model # LBP-742D

1. Battery voltage: 7.4 V DC
2. Battery capacity: 1900 mAh
3. Output power: max. 15 W

NOTE: Do not connect the rechargeable battery to other batteries or charging devices which are not part of the delivery.

4. Battery charge: During the charging process, the battery LED is lit in red. It will change to green as soon as the battery is fully charged. The battery can be recharged about 500 times.

C. Product description

1. Voltage: 7.4 V
2. After setting the device to "High" (LED is lit in red), the battery can be used for 3 hours. After setting the device to "Low" (LED is lit in green), the battery can be used for 6 hours.

The functionality of this item has been tested before the delivery.

Please do not rely on simply assessing the temperature with your hand as the hands are not as sensitive to temperatures as the skin surface in other parts of the body. Therefore a temperature that feels lukewarm on your hands may feel unpleasantly hot on the skin elsewhere.

As infrared radiation primarily warms up the surfaces while having only a secondary warming effect on the air, the air temperature may be situated well below the temperature felt by a person. In contrast to conventional kinds of heat therapy which only warm up the skin surface, infrared warmth reaches well into the muscles and bones as well. The heat efficiency of infrared radiation enables the user to apply this therapy at an all in all lower temperature level – which means it is ideally suited for heat-sensitive persons.

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